



MEDICAL-X
ADVANCING SIMULATION



IV ARM

Task Trainer

BROCHURE

IV Arm

- Advanced arm model for venipuncture practice
- Life-like touch and feel
- Realistic and faithful anatomy
- Durable material and leakage-free
- Easy to assemble and transport

Meet the IV Arm

The Medical-X's IV Arm is the most realistic and faithful arm simulator for practicing peripheral intravenous cannulation and blood sampling.

With an accurate venous system of the upper limb and several underlying palpable veins, the IV Arm also has a soft silicone coat that reproduces the feeling of touching and puncturing the human skin, increasing the realism of your training.

An add-on deltoid is available to perform intramuscular and subcutaneous injections, completing the training session.

Multiple sites for venipuncture




- Basilic Vein
- Cephalic Vein
- Antecubital Vein
- Dorsal Metacarpal Veins

Life-size arm model with supporting stand and closed blood flow system




Realistic blood flashback




IV cannulation



Deltoid Add-on for IM and SC injections



Replaceable veins



Specification list

- Life-size arm model for venipuncture training
- Realistic soft and flexible silicone skin
- Faithful anatomy of the upper limb venous system
- Palpable veins on the arm and hand with realistic tactile feeling
- Multiple sites for practicing: Basilic Vein, Cephalic Vein, Antecubital Vein, and Dorsal Metacarpal Veins
- Easy and quick assembling
- A closed blood flow system keeps the veins' pressure
- Realistic puncture resistance
- Realistic fluid flash-back into the injection needle
- Can be used with a vacuum blood collection system
- Needle catheter easily slides during vein cannulation
- Durable and resistant skin
- Durable blood veins
- Leak-resistant veins
- Replaceable skin
- Veins are easily accessible and easy to replace
- Easy to clean and disinfect
- Lightweight unit and easy to carry-on
- Includes a carrying bag
- A Deltoid add-on can be purchased separately and attached to the IV Arm for intramuscular and subcutaneous injections training

General Warning:

The **IV Arm** and its components contain silicone and latex. Do not use this product if you have any known allergy to these materials.

Discontinue use if skin irritation occurs.

The **IV Arm** contains delicate parts. Handle with care to prevent damage.

It is recommended to use needles/catheters smaller than 20G to ensure the product's longevity.



MEDICAL-X
ADVANCING SIMULATION

Scan QR code
to discover more



medical-x.com | info@medical-x.com | +31 10 34 15 927

Schiedamsedijk 77 | 3011 EM, Rotterdam